



Technical buffer solutions

Technical buffer solutions

Identified by colour-coded labels. Tolerance ± 0.03 pH.



Code	Description	Packaging
LLG09040722	Technical buffer solutions, pHvalue 4.01 at 25 °C, 1000 ml	1 pz.
LLG09040723	Technical buffer solutions, pHvalue 7.00 at 25 °C, 1000 ml	1 pz.
LLG09040725	Technical buffer solutions, pHvalue 10.01* at 25 °C, 1000 ml	1 pz.
LLG09040727	Technical buffer solutions, pHvalue 4.01 at 25 °C, 250 ml	1 pz.
LLG09040728	Technical buffer solutions, pHvalue 7.00 at 25 °C, 250 ml	1 pz.
LLG09040734	Technical buffer solutions, pHvalue 4.01 at 25 °C, 50 ml	1 pz.
LLG09040735	Technical buffer solutions, pHvalue 7.00 at 25 °C, 50 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

