




Positioning lids for Sonorex ultrasonic baths

Bandelin

Positioning lids for Sonorex ultrasonic baths

Stainless steel.

For supporting inset beakers.



Code	Description	Packaging
LLG09877617	Positioning lids for Sonorex ultrasonic baths, Type ES 4, For bath type RK 1028/H/CH, DT 1028 / H/CH, RK 1050 / DT 1050 CH, DL 1028 H	1 pz.
LLG09877613	Positioning lids for Sonorex ultrasonic baths, Type DE 100, For bath type RK/DT 100 H / 102 H / 103 H, DL 102 H	1 pz.
LLG09877616	Positioning lids for Sonorex ultrasonic baths, Type DE 6, For bath type RK/DT 106	1 pz.
LLG09877612	Positioning lids for Sonorex ultrasonic baths, Type DE 156, For bath type RK/DT 156 BH / 156 BH-RC, DL 156 BH	1 pz.
LLG09877615	Positioning lids for Sonorex ultrasonic baths, Type DE 255, For bath type RK/DT 255 H / 255 H-RC, DL 255 H	1 pz.
LLG09877621	Positioning lids for Sonorex ultrasonic baths, Type DE 08, For bath type RK/DT 31 / H	1 pz.
LLG09877620	Positioning lids for Sonorex ultrasonic baths, Type DE 510, For bath type RK/DT 510 / 510 H / 512 H, DL 510 H	1 pz.
LLG09877614	Positioning lids for Sonorex ultrasonic baths, Type DE 514, For bath type RK/DT 514 H / 514 BH / 514 BH-RC, DL 514 BH	1 pz.
LLG09877611	Positioning lids for Sonorex ultrasonic baths, Type DE 52, For bath type RK/DT 52 / 52H / 52H-RC	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

