


**Positioning lids for Sonorex ultrasonic baths****Bandelin****Positioning lids for Sonorex ultrasonic baths**

Stainless steel.

For supporting inset beakers.



Code	Description	Packaging
LLG09877617	Positioning lids for Sonorex ultrasonic baths, Type ES 4, For bath type RK 1028/H/CH, DT 1028 / H/CH, RK 1050 / DT 1050 CH, DL 1028 H	1 pz.
LLG09877613	Positioning lids for Sonorex ultrasonic baths, Type DE 100, For bath type RK/DT 100 H / 102 H / 103 H, DL 102 H	1 pz.
LLG09877616	Positioning lids for Sonorex ultrasonic baths, Type DE 6, For bath type RK/DT 106	1 pz.
LLG09877612	Positioning lids for Sonorex ultrasonic baths, Type DE 156, For bath type RK/DT 156 BH / 156 BH-RC, DL 156 BH	1 pz.
LLG09877615	Positioning lids for Sonorex ultrasonic baths, Type DE 255, For bath type RK/DT 255 H / 255 H-RC, DL 255 H	1 pz.
LLG09877621	Positioning lids for Sonorex ultrasonic baths, Type DE 08, For bath type RK/DT 31 / H	1 pz.
LLG09877620	Positioning lids for Sonorex ultrasonic baths, Type DE 510, For bath type RK/DT 510 / 510 H / 512 H, DL 510 H	1 pz.
LLG09877614	Positioning lids for Sonorex ultrasonic baths, Type DE 514, For bath type RK/DT 514 H / 514 BH / 514 BH-RC, DL 514 BH	1 pz.
LLG09877611	Positioning lids for Sonorex ultrasonic baths, Type DE 52, For bath type RK/DT 52 / 52H / 52H-RC	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

