

**Vacuum Filtration Equipment GV 050 Series****Whatman****Vacuum Filtration Equipment GV 050 Series**

- Funnel Volume 250 ml
- to be used with Ø 47/50 mm membranes and Ø 40 mm pre-filter
- silicone cap and venting nozzle included in scope of supply
- GV 050A-Types equipped with rapid closure clamp



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09603500 | Vacuum filtration equipment GV 050/050m | 1 pz. |
| LLG09603502 | Vacuum filtration equipment GV 050/250m | 1 pz. |
| LLG09603501 | Vacuum filtration equipment GV 050/150m | 1 pz. |
| LLG09603503 | Vacuum filtration equipment GV 050/350m | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

