

**Vacuum filtration Equipment Grade MV 050/0****Whatman****Vacuum filtration Equipment Grade MV 050/0**

Stainless steel construction, chemically resistant to most aqueous and organic solutions.

- Funnel volume 500 ml
- For filters with 47 mm /50 mm diameter and 40 mm diameter prefilter
- Self sealing, ideal for quick, repeated filtrations



Code	Description	Packaging
LLG09603510	Apparecchio per filtrazioni sottovuoto	1 pz.
LLG09603520	Vacuum filtration units,metal,for 47/50	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

