

**Narrow-mouth bottles with screw thread****VITLAB****Narrow-mouth bottles with screw thread**

PFA. Transparent, with tightly sealing screw closure (saw-tooth thread). High temperature resistance from -200 to +250°C. Autoclavable at 121°C (2 bar) according to DIN EN 285.



Code	Description	Packaging
LLG09102420	Narrow-mouth bottles with screw thread, 100 ml, Diam. 45 mm, Height 120 mm, Thread 28 S	1 pz.
LLG09102426	Narrow-mouth bottles with screw thread, 1000 ml, Diam. 96 mm, Height 240 mm, Thread 28 S	1 pz.
LLG09102422	Narrow-mouth bottles with screw thread, 250 ml, Diam. 61 mm, Height 160 mm, Thread 28 S	1 pz.
LLG09073740	Narrow-mouth bottles with screw thread, 50 ml, Diam. 37 mm, Height 86 mm, Thread 28 S	1 pz.
LLG09102424	Narrow-mouth bottles with screw thread, 500 ml, Diam. 76 mm, Height 190 mm, Thread 28 S	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

