

**Sample tubes, PFA****VITLAB****Sample tubes, PFA**

Tube for sample preparation and for use in autosampler racks. With or without individually calibrated ring mark at 10 ml with GL 25 screw cap made from PFA or PE stoppers.



Code	Description	Packaging
LLG09315345	Sample tubes, PFA, with ring mark and screw cap, 15 ml, Diam. 22 mm, Height 110 mm	1 pz.
LLG09315344	Sample tubes, PFA, with ring mark and stopper, 12 ml, Diam. 16 mm, Height 110 mm	1 pz.
LLG09315346	Sample tubes, PFA, without ring mark and screw cap, 15 ml, Diam. 22 mm, Height 110 mm	1 pz.
LLG09315343	Sample tubes, PFA, without ring mark and stopper, 12 ml, Diam. 16 mm, Height 110 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

