

**Gas washing bottles, PFA****VITLAB****Gas washing bottles, PFA**

With cap with S40 screw thread and srewable frit made of PTFE. Approx. pore size 3 μm for optimal formation of gas bubbles into the liquid. Because all parts are made of high-quality fluorocarbon resins, a broad range of applications are possible. The threaded inlet and outlet connections accept 4 x 6mm i.d. x o.d. tubing , (1000 ml version: 5 x 8mm i.d. x o.d. tubing).



Code	Description	Packaging
LLG09110324	Gas washing bottles, PFA, 1000 ml, Diam. 96 mm, Height 240 mm	1 pz.
LLG09110320	Gas washing bottles, PFA, 250 ml, Diam. 61 mm, Height 160 mm	1 pz.
LLG09110322	Gas washing bottles, PFA, 500 ml, Diam. 76 mm, Height 190 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

