



Buffer solutions, sterile

SI Analytics

Buffer solutions, sterile

According to DIN 19 266. Values stated at 25°C. To be delivered in DURAN® glass bottles. Hot steam sterilized. Accuracy ± 0.012 pH. Comes with individual batch number and manufacturer certificate.



Code	Description	Packaging
LLG09041007	Buffer solutions, sterile, pHvalue 6.87 at 25 °C, 250 ml	1 pz.
LLG09041004	Buffer solutions, sterile, pHvalue 4.01 at 25 °C, 250 ml	1 pz.
LLG09041009	Buffer solutions, sterile, pHvalue 9.18 at 25 °C, 250 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

