


**Buffer solutions, sterile****SI Analytics****Buffer solutions, sterile**

According to DIN 19 266. Values stated at 25°C. To be delivered in DURAN® glass bottles. Hot steam sterilized. Accuracy  $\pm 0.012\text{pH}$ . Comes with individual batch number and manufacturer certificate.



Code	Description	Packaging
LLG09041007	Buffer solutions, sterile, pHvalue 6.87 at 25 °C, 250 ml	1 pz.
LLG09041004	Buffer solutions, sterile, pHvalue 4.01 at 25 °C, 250 ml	1 pz.
LLG09041009	Buffer solutions, sterile, pHvalue 9.18 at 25 °C, 250 ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

