



## Single channel microliter pipettes Eppendorf Research plus, variable

### Single channel microliter pipettes Eppendorf Research® plus, variable

Available with fixed or adjustable volume. The ultra light Eppendorf Research plus pipettor meets the highest requirements in precision and accuracy.

Ergonomic:

- Color-coded control buttons for volume detection
- Four-digit volume display with magnifying display
- A few turns to set the desired volume
- Very low operating forces (control knob, ejector)
- Spring-loaded tip cone for optimum tip fit with minimal attachment force

Flexible:

- Selection between single-channel, multi-channel and fixed volume pipettes in different sizes
- Adjustment possibilities for liquids of different density than water

Robust:

- High resistance to chemicals and mechanical stress
- Fully autoclavable

-Eppendorf PerfectPistonä: Ultra light piston made of Fortron® (except 2.5µl, 10µl, 5ml and 10ml pipettors)

-quick connection clip: easily removable lower part (not in 5ml and 10ml pipettes)

-spring loaded tip cone: Improved ergonomics, tight fit to the tip (except 5ml and 10ml pipettors)





**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

