



Pipette tips Qualitix, natural, thin, 10 µl

Pipette tips Qualitix®, natural, thin, 10 µl

The Socorex Qualitix® are exclusively made from high grade, virgin, metal-free polypropylene. Precision injection moulding care for air-tight fit. Stringent quality procedures ensure lot by lot consistency.

- Wide selection, superior quality
- High density filter prevents aerosol contamination
- DNAse, RNAse, ATP and pyrogen-free tips
- Compatibility with various pipette brands



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09920715 | PTFE collars PTFE 29for RG 2, SR 2 and | 12 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

