



Incubators, BD, BF series - Classic.Line

Incubators, BD, BF series - Classic.Line


BINDER incubators stand for optimal incubation and for results that are reproducible in every routine test, even under high batch throughputs in long-term operation.

Serie BD: Incubators with gravity convection.

Serie BF: Incubators with forced convection.

Equipment:

- Temperature range: from 5°C above room temperature to 100°C
- Independent temperature safety device class 3.1 (DIN 12880) with optical and audible temperature alarm
- RS 422 interface for use with APT-COMä DataControlSystem communication software
- Inner door made of tempered safety glass
- Power supply 230 V, 50/60 Hz



Code	Description	Packaging
LLG09883599	Incubators BD 23 21l, with RS 422	1 pz.
LLG09883603	Incubator, Model BD 400, upto 99,9°C, cap. 400 ltrs	1 pz.
LLG09883572	Cooled incubator KB 720 standard, 698 l, 200...240 V 1~ 50/60 Hz, with RD3.1	1 pz.
LLG09883570	Cooled incubator KB 240 standard, 247 l, 200...240 V 1~ 50/60 Hz, with RD3.1	1 pz.
LLG09883571	Cooled incubator KB 400 standard, 400 l, 200...240 V 1~ 50/60 Hz, with RD3.1	1 pz.
LLG09883528	Cooled incubator KB 23 Standard, 23 Ltr. 230V. -50/60 Hz-, with RD3.1	1 pz.
LLG09883569	Cooled incubator KB 115 standard, 115 l, 230 V 1~ 50 Hz, with RD3.1	1 pz.
LLG06203549	Incubator BF 400 with forced convection, 400 l 230 V, 50/60Hz	1 pz.
LLG09883568	Cooled incubator KB 053 standard, 53 l, 230 V 1~ 50 Hz, with RD3.1	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

