



Tubing Tygon 3350

Tubing Tygon® 3350

Silicone-platinum tubing, Shore A 50. For sterile filling systems, chemistry and blood analyzers and liquid chromatography. Biocompatibility for sensitive applications. To FDA and USP Class VI, in accordance with ISO 10993.

- Temperature range -81 to 200°C
- Smooth inner surface
- High flexibility
- Ozone and UV resistant



Code	Description	Packaging
LLG09208383	Tubing Tygon® 3350, Int. diam. 6.4 mm, Ext. diam. 11.2 mm, Wall thickness 2.4 mm, Shore 50, Temp. range -81 ... 200 °C	15 m
LLG09208385	Tubing Tygon® 3350, Int. diam. 9.5 mm, Ext. diam. 12.7 mm, Wall thickness 1.6 mm, Shore 50, Temp. range -81 ... 200 °C	15 pz
LLG09208384	Tubing Tygon® 3350, Int. diam. 8 mm, Ext. diam. 12.8 mm, Wall thickness 2.4 mm, Shore 50, Temp. range -81 ... 200 °C	15 m
LLG09208386	Tubing Tygon® 3350, Int. diam. 12.7 mm, Ext. diam. 15.9 mm, Wall thickness 1.6 mm, Shore 50, Temp. range -81 ... 200 °C	15 pz
LLG09208387	Tubing Tygon® 3350, Int. diam. 12.7 mm, Ext. diam. 17.5 mm, Wall thickness 2.4 mm, Shore 50, Temp. range -81 ... 200 °C	15 pz
LLG09208380	Tubing Tygon® 3350, Int. diam. 3.2 mm, Ext. diam. 6.4 mm, Wall thickness 1.6 mm, Shore 50, Temp. range -81 ... 200 °C	15 pz
LLG09208381	Tubing Tygon® 3350, Int. diam. 4.8 mm, Ext. diam. 8 mm, Wall thickness 1.6 mm, Shore 50, Temp. range -81 ... 200 °C	15 pz
LLG09208382	Tubing Tygon® 3350, Int. diam. 6.4 mm, Ext. diam. 9.6 mm, Wall thickness 1.6 mm, Shore 50, Temp. range -81 ... 200 °C	15 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
 Reagents operates with
 a Certified Quality
 Management System

