



Beakers, Borosilicate glass 3.3, tall form

Beakers, Borosilicate glass 3.3, tall form

- DIN 12331
- ISO 3819
- With approximate graduations and spout



Code	Description	Packaging
LLG06270117	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 2000 ml	10 pz
LLG09013924	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 25 ml *	10 pz.
LLG09013925	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 50 ml	10 pz.
LLG09013920	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 100 ml	10 pz.
LLG09013921	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 150 ml	10 pz.
LLG09013922	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 250 ml	10 pz.
LLG09013923	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 400 ml	10 pz.
LLG09013930	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 600 ml	10 pz.
LLG09013931	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 800 ml	10 pz.
LLG09013932	Beakers, Borosilicate glass 3.3, tall form, Nominal capacity 1000 ml	10 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

