

**Stereomicroscope, ST-36C****Motic****Stereomicroscope, ST-36C**

- Binocular head, 45° inclined, 360° rotating
- Widefield eyepieces WF10X/20mm
- Objective turret magnification changer (2X, 4X)
- Coarse focusing system with tension adjustment
- Base stand with pole
- Incident and transmitted rechargeable LED illumination (ST36C-6LED)
- Incident and transmitted illumination 12V/10W (ST36C-2L00)
- Power supply 100V-240V (CE), charger and batteries included (ST36C-6LED)
- Main supply 220V-240V (CE) (ST36C-2L00)
- Black/white and frosted glass stage plates Ø 80mm, dust cover



Code	Description	Packaging
LLG09727007	Stereomicroscope, ST-36C, Type ST-36C-2LOO,, Lighting Tungsten	1 pz
LLG09727008	Stereomicroscope, ST-36C, Type ST-36C-2LOO UK,, Lighting Tungsten	1 pz
LLG09727009	Stereomicroscope, ST-36C, Type ST-36C-6LED, battery, Lighting LED	1 pz
LLG09727010	Stereomicroscope, ST-36C, Type ST-36C-6LED UK, battery, Lighting LED	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

