


**Electronic multi-size heating mantles, V-shaped, EMV series****Electrothermal****Electronic multi-size heating mantles, V-shaped, EMV series**

Solid, with mesh stainless steel liner, which protects against fluid spills, with central aperture to accommodate round-bottom and pear-shaped flask and 60° funnels. Polypropylene outer case is resilient and chemically-resistant. Highly efficient in heating up to a maximum element temperature of 450°C. Even at full power output, the exterior remains "cool-to-touch" due to good heating element insulation. All heating mantles have support clamp for rods of up to 12mm diameter, and are double fused, with earth (ground) screen to protect the user from electric shocks. Available as 220-240 V 50/60 Hz models; spare heater cartridges are available on request.

For round bottom flasks, pear-shaped flasks and 60° funnels (Mantle is open at the bottom).



Code	Description	Packaging
LLG09643072	Electronic multi-size heating mantles, V-shaped, EMV series, For flasks 10 to 50 ml, For funnel diam. 50 to 100 mm	1 pz.
LLG09643074	Electronic multi-size heating mantles, V-shaped, EMV series, For flasks 500 to 1000 ml, For funnel diam. 100 to 200 mm	1 pz.
LLG09643073	Electronic multi-size heating mantles, V-shaped, EMV series, For flasks 100 to 250 ml, For funnel diam. 75 to 100 mm	1 pz.
LLG09643075	Electronic multi-size heating mantles, V-shaped, EMV series, For flasks 2000 to 5000 * ml, For funnel diam. 200 to 300 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

