

**Electronic stirrer mantles, EMA series****Electrothermal****Electronic stirrer mantles, EMA series**

Built-in adjustable heater control up to an element temperature of 450°C and magnetic stirring control up to approx. 520rpm with auto-reverse facility and automatic stirrer bar trap; the stirrer unit has an independent power supply. Polypropylene outer case is resilient and chemically-resistant. Even at full power output, the exterior remains "cool-to-touch" due to good heating element insulation. All heating mantles have support clamp for rods of up to 12mm diameter, and are double fused, with earth (ground) screen to protect the user from electric shocks. Available as 220-240 V 50/60 Hz models; spare heater cartridges are available on request.



Code	Description	Packaging
LLG09643124	Electronic stirrer mantles, EMA series, For flasks 250 ml, Width 175 mm, Depth 260 mm, Height 133 mm, Power 170 W	1 pz.
LLG09643122	Electronic stirrer mantles, EMA series, For flasks 50 ml, Width 175 mm, Depth 260 mm, Height 133 mm, Power 80 W	1 pz.
LLG09643123	Electronic stirrer mantles, EMA series, For flasks 100 ml, Width 175 mm, Depth 260 mm, Height 133 mm, Power 80 W	1 pz.
LLG09643125	Electronic stirrer mantles, EMA series, For flasks 500 ml, Width 238 mm, Depth 310 mm, Height 157 mm, Power 220 W	1 pz.
LLG09643126	Electronic stirrer mantles, EMA series, For flasks 1000 ml, Width 238 mm, Depth 310 mm, Height 157 mm, Power 320 W	1 pz.
LLG09643127	Electronic stirrer mantles, EMA series, For flasks 2000 ml, Width 350 mm, Depth 400 mm, Height 197 mm, Power 520 W	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

