

**Hand spectromscopes****A.KRUSS Optronic****Hand spectromscopes**

Handy and precise. Angle dispersion C -F: 7°. Linear dispersion: 60 mm (direct vision prism of Amici type)

Ideal for demonstrations in schools and for:

- Viewing / observation of emission and absorption spectra
- Control of the wavelength spectra of gas emission lamps



Code	Description	Packaging
LLG09841501	Hand spectrometers, Type 1501, Labeling with adjustable split	1 pz.
LLG09841504	Hand spectrometers, Type 1504, Labeling with adjustable split and wavelength scale (400 to 750nm)	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

