

**Desiccator stopcocks, borosilicate glass 3.3****Lenz****Desiccator stopcocks, borosilicate glass 3.3**

With PTFE spindle. For NOVUS-type desiccators.

Cone: NS 24/29, bore: NW 2,5mm,, hose connection: outer dia. 8mm.



Code	Description	Packaging
LLG09042505	Desiccator stopcocks, borosilicate glass 3.3, for lid tubulation, Length 124 mm, NS 24/29	1 pz.
LLG09042500	Desiccator stopcocks, borosilicate glass 3.3, for side tubulation, Length 108 mm, NS 24/29	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

