

**Magnetic stirrer, Cimarec i Maxi Direct****Thermo Scientific****Magnetic stirrer, Cimarec i Maxi Direct**

With non-mechanical, electronic induction drive and dust / watertight, easy-clean stainless steel housing. Constant speed irrespective of load even if viscosity changes. With large footprint and a strong and more penetrating magnetic field. Not suitable for use in incubators. With mains adapter.

**Specifications**

**Max. volume stirred (H<sub>2</sub>O): 5 l**

**Motor rating input: 20 W**

**Stirring positions: 1**

**Speed range: 80 $\frac{1}{4}$ 2000 min<sup>-1</sup>**

**Permissible ambient temperature: -10 $\frac{1}{4}$ 40 °C**

**Permissible humidity: 95 %**

**Display: digital**

**Dimensions (W x D x H): 240 x 215 x 35 mm**

**Weight: 1,4 kg**

**IP code: IP 64**

**Power supply: 115/230 V, 50/60 Hz**



| Code        | Description                                                                                    | Packaging |
|-------------|------------------------------------------------------------------------------------------------|-----------|
| LLG09720618 | Magnetic stirrer MAXI DIRECT,with power supply 180 x 180 x 35 mm,stirring volume:250 - 4000 ml | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

