

**Multi-position magnetic stirrers, Cimarec i Poly 15****Thermo Scientific****Multi-position magnetic stirrers, Cimarec i Poly 15**

Without heating. With non-mechanical, electronic induction drive. With 15 stirring spaces, for e.g. 6 x 1L erlenmeyer flasks or 15 x 250ml beakers. For mass screening and biological applications, even in humidity incubators. With built-in controller, digital speed indicator and soft-start drive. With mains adaptor.

Specifications

Max. stirring volume (H₂O): 3 l

Motor rating input: 10 W

Stirring positions: 15

Speed range: 130 ... 990 min⁻¹

Permissible ambient temperature: -10 ... 40 °C

Permissible humidity: 95 %

Display: digital

Dimensions (W x D x H): 240 x 420 x 35 mm

Weight: 6 kg

IP code: IP 32

Power supply: 115/230 V, 50/60 Hz



Code	Description	Packaging
LLG09720650	Magnetic stirrer POLY 15 grey, with 15 stirring positions speed: 100-980 rpm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

