

**Multi-position magnetic stirrers, Cimarec i Poly 15****Thermo Scientific****Multi-position magnetic stirrers, Cimarec i Poly 15**

Without heating. With non-mechanical, electronic induction drive. With 15 stirring spaces, for e.g. 6 x 1L erlenmeyer flasks or 15 x 250ml beakers. For mass screening and biological applications, even in humidity incubators. With built-in controller, digital speed indicator and soft-start drive. With mains adaptor.

**Specifications**

**Max. stirring volume (H<sub>2</sub>O): 3 l**

**Motor rating input: 10 W**

**Stirring positions: 15**

**Speed range: 130 ... 990 min<sup>-1</sup>**

**Permissible ambient temperature: -10 ... 40 °C**

**Permissible humidity: 95 %**

**Display: digital**

**Dimensions (W x D x H): 240 x 420 x 35 mm**

**Weight: 6 kg**

**IP code: IP 32**

**Power supply: 115/230 V, 50/60 Hz**



| Code        | Description                                                                  | Packaging |
|-------------|------------------------------------------------------------------------------|-----------|
| LLG09720650 | Magnetic stirrer POLY 15 grey, with 15 stirring positions speed: 100-980 rpm | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

