

**LLG single channel microliter pipettes, variable****LLG Labware****LLG single channel microliter pipettes, variable****Available with fixed or adjustable volume.****Ergonomic shape for a very comfortable handling**

The rounded shape and the unique soft grip of the LLG pipette ensures a safe and comfortable handling. A minimal effort is required for pipetting so that the risk of repetitive strain injuries (RSI) is reduced.

Maximum Precision and Accuracy

Thermal insulation of the internal components, using the innovative soft grip, improves the accuracy.

Universal Blow-out System

The two-step plunger operation allows reverse pipetting. This also facilitates last drop dispensing by soft plunger movements. High accuracy and precision are guaranteed.

Colour-coded push-button for quick operation

The colour code enables a quick choice of the right pipette depending of the volume to be pipetted (see corresponding column in the table below).

Easy Maintenance and Adjustment

The volume setting is easy and precise, using the push-button, even when wearing gloves. Maintenance and adjustment take only a few minutes with just one tool. (Included in delivery)

Contamination-free Pipetting

The pipette is fully autoclavable, ideal for applications where sterility is important.

Satisfaction Guaranteed

Each pipette is individually tested and calibrated twice, conforming to EN ISO 8655 and DIN 12600 Standards. A certificate is included.



Made in Germany, CE certified.

Code	Description	Packaging
LLG09280004	LLG single channel microliter pipettes, variable, 20 - 200 μ l, Increments 1.0 μ l	1 pz
LLG09280005	LLG single channel microliter pipettes, variable, 100 - 1000 μ l, Increments 5.0 μ l	1 pz
LLG09280007	LLG single channel microliter pipettes, variable, 1000 - 10000* μ l, Increments 100.0 μ l	1 pz
LLG09280003	LLG single channel microliter pipettes, variable, 10 - 100 μ l, Increments 1.0 μ l	1 pz



Code	Description	Packaging
LLG09280002	LLG single channel microliter pipettes, variable, 2 - 20 µl, Increments 0.5 µl	1 pz
LLG09280001	LLG single channel microliter pipettes, variable, 0.5 - 10 µl, Increments 0.1 µl	1 pz
LLG09280000	LLG single channel microliter pipettes, variable, 0.1 - 2.5 µl, Increments 0.05 µl	1 pz
LLG09280006	LLG single channel microliter pipettes, variable, 500 - 5000 µl, Increments 50.0 µl	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

