

**Disposable liquid samplers, HDPE****Bürkle****Disposable liquid samplers, HDPE**

In numerous situations, rapid and contamination-free samples have to be taken. DispoSamplers are perfect for this application. They are cleanroom-manufactured, individually-packed and always ready-to-use. Avoid expensive decontamination or cleaning of soiled samplers. Gamma-sterilised versions are also available.

ViscoDispo: For high viscosity-media up to 100000 mPas.

LiquiDispo: For liquids and slightly viscous media.

DispoPipette: For taking all-layer and target-point samples from low viscosity-media.

DispoTube: For taking all-layer and target-point samples from liquids.



Code	Description	Packaging
LLG06282447	Disposable liquid samplers, HDPE, Type DispoTube, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG06282448	Disposable liquid samplers, HDPE, Type DispoTube, sterile, 100 ml, Length 500 mm	20 pz
LLG06282451	Disposable liquid samplers, HDPE, Type DispoPipette, non-sterile, 100 ml, Length 500 mm	20 pz.
LLG07971434	Disposable liquid samplers, HDPE, Type DispoPipette, sterile, 100 ml, Length 500 mm	20 pz
LLG09303721	Disposable liquid samplers, HDPE, Type LiquiDispo, sterile, 100 ml, Length 500 mm	20 pz.
LLG09303733	Disposable liquid samplers, HDPE, Type ViscoDispo, sterile, 100 ml, Length 500 mm	20 pz
LLG06281401	Disposable liquid samplers, HDPE, Type DispoPipette, sterile, 200 ml, Length 1000 mm	20 pz
LLG06281402	Disposable liquid samplers, HDPE, Type DispoPipette, non-sterile, 200 ml, Length 1000 mm	20 pz.
LLG06282450	Disposable liquid samplers, HDPE, Type DispoTube, sterile, 200 ml, Length 1000 mm	20 pz
LLG09303723	Disposable liquid samplers, HDPE, Type LiquiDispo, sterile, 200 ml, Length 1000 mm	20 pz
LLG09303735	Disposable liquid samplers, HDPE, Type ViscoDispo, sterile, 200 ml, Length 1000 mm	20 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

