

**LLG-Pipette controller easy****LLG Labware****LLG-Pipette controller easy**

The LLG-Pipette controller easy enables fast and exact pipetting with glass and plastic pipettes.

- Ergonomic and lightweight construction for constant ease of use over prolonged usage
- Unique nose cone angle for enhanced flexibility in operation e.g. in fume hoods
- Eco-friendly rechargeable NiMH batteries, offering up to 8 hours of continuous working
- Easily replaceable batteries
- Autoclavable tip cone for avoiding cross-contamination

**Scope of supply:** Pipette controller, manual, charger, wall-mount, desktop stand



Code	Description	Packaging
LLG06263683	LLG-pipette controller easy, EU-plug	1 pz
LLG06263684	LLG-pipette controller easy, UK-plug	1 pc



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

