

**Homogenisers Potter-Elvehjem, with glass pestle****DWK Life Sciences****Homogenisers Potter-Elvehjem, with glass pestle**

All-glass. Pestle is suitable for motor drive. The space between the pestle and mortar is between 0.1mm to 0.15mm. Spare components are available and are completely interchangeable.



Code	Description	Packaging
LLG09651677	Homogenisers Potter-Elvehjem, with glass pestle, 17 ml, Size 23, Pestle L x o.d. 265 x 10 mm, Mortar L x o.d. 175 x 25 mm	1 pz
LLG09651676	Homogenisers Potter-Elvehjem, with glass pestle, 8 ml, Size 22, Pestle L x o.d. 215 x 8 mm, Mortar L x o.d. 150 x 18 mm	1 pz
LLG09651678	Homogenisers Potter-Elvehjem, with glass pestle, 45 ml, Size 24, Pestle L x o.d. 315 x 10 mm, Mortar L x o.d. 220 x 32 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

