

**Precision balances PES / PEJ****Kern & Sohn****Precision balances PES / PEJ**

Robust laboratory and industrial precision scale for heavy items.

- Only PEJ: Automatic internal adjustment
- Only PES: Adjusting program CAL for quick setting of the balance accuracy, with an external test weight (must be purchased separately)
- Metal housing: robust and sturdy
- Weighing with tolerance range (checkweighing): Input of 2 upper and 2 lower limit values through 4 arrow keys. An audible and visual signal assists with the portioning, dispensing or grading
- Draught shield standard on all models with readout = 0,001 g, weighing space (W x D x H): 170 x 150 x 100 mm. Removable metal cover with pipette opening
- Underfloor weighing: load support with hook on the underside of the balance (only for models with a weighing plate size of 140x120 mm and 200x200 mm)
- Fluorescent display, bright with high contrast, digit height 13 mm
- Permissible ambient temperature: 10 °C / 30 °C



Code	Description	Packaging
LLG07651840	Precision balances PES 620-3M, 620 g	1 pz.
LLG07670597	Precision balances PEJ 620-3M*, 620 g	1 pz.
LLG06237670	Precision balances PES 31000-1M, 31000 g	1 pz.
LLG07920074	Precision balances PES 15000-1M, 15000 g	1 pz.
LLG06225841	Precision balances PES 2200-2M, 2200 g	1 pz.
LLG06229872	Precision balances PES 4200-2M, 4200 g	1 pz.
LLG06235603	Precision balances PEJ 4200-2M*, 4200 g	1 pz.
LLG06236499	Precision balances PES 6200-2M, 6200 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

