

**Thermostatic baths with polycarbonate tanks****Huber****Thermostatic baths with polycarbonate tanks**

The transparent polycarbonate baths are suitable for use up to 100°C. An immersion thermostat is mounted on the bath bridge for all models. With a pump adapter, this combination can also be used with external, closed applications. Models with the Pilot ONE controller have a variable speed pressure/suction pump and are therefore also suitable for external open applications.

Temperature max.: 100 °C

Temp. Stability at 70 °C to DIN12876

KISS E:

CC-E: 0.05K

0.02K

Safety class: FL, III

Bath depth: 150mm



Code	Description	Packaging
LLG06272334	Thermostatic baths with polycarbonate tanks, Type KISS-106A, Bath capacity 6 L, Pressure pump max. 14 / 0.25 L/min bar	1 pz.
LLG09859203	Thermostatic baths with polycarbonate tanks, Type CC-106A, Bath capacity 6 L, Pressure pump max. 27 / 0.7 L/min bar	1 pz.
LLG06272335	Thermostatic baths with polycarbonate tanks, Type KISS-108A, Bath capacity 8 L, Pressure pump max. 14 / 0.25 L/min bar	1 pz.
LLG09859204	Thermostatic baths with polycarbonate tanks, Type CC-108A, Bath capacity 8 L, Pressure pump max. 27 / 0.7 L/min bar	1 pz.
LLG06272337	Thermostatic baths with polycarbonate tanks, Type KISS-112A, Bath capacity 12 L, Pressure pump max. 14 / 0.25 L/min bar	1 pz.
LLG09859206	Thermostatic baths with polycarbonate tanks, Type CC-112A, Bath capacity 12 L, Pressure pump max. 27 / 0.7 L/min bar	1 pz.
LLG06272338	Thermostatic baths with polycarbonate tanks, Type KISS-118A, Bath capacity 18 L, Pressure pump max. 14 / 0.25 L/min bar	1 pz.
LLG09859207	Thermostatic baths with polycarbonate tanks, Type CC-118A, Bath capacity 18 L, Pressure pump max. 27 / 0.7 L/min bar	1 pz.
LLG06272336	Thermostatic baths with polycarbonate tanks, Type KISS-110A, Bath capacity 10 L, Pressure pump max. 14 / 0.25 L/min bar	1 pz.
LLG09859205	Thermostatic baths with polycarbonate tanks, Type CC-110A, Bath capacity 10 L, Pressure pump max. 27 / 0.7 L/min bar	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

