



Centrifugal concentrators, Vivaspin 2

Centrifugal concentrators, Vivaspin 2

Membrane: polyethersulphone. This 2 ml volume unit combines the rapidity of the classic Vivaspin products with reduced internal surface and membrane areas for superior recoveries from very dilute solutions. With a choice of suitable membranes (PES, CTA, Hydrosart) Vivaspin 2 provides the highest flexibility for process optimisation. Vivaspin 2 guarantees maximum sample recovery either by reverse centrifuging into the concentrate recovery cap or by directly pipetting the concentrate from the deadstop pocket built into the bottom of the concentrator.

Membrane surface: 1.2 cm²

Max. rcf. : 8000 x g

Hold-up volume: 10 µl



Code	Description	Packaging
LLG09049861	Centrifugal concentrators, Vivaspin 2, MWCO 3000 da	25 pz.
LLG09049862	Centrifugal concentrators, Vivaspin 2, MWCO 3000 da	100 pz.
LLG09049760	Centrifugal concentrators, Vivaspin 2, MWCO 5000 da	25 pz.
LLG09049777	Centrifugal concentrators, Vivaspin 2, MWCO 5000 da	100 pz.
LLG09049761	Centrifugal concentrators, Vivaspin 2, MWCO 10000 da	25 pz.
LLG09049778	Centrifugal concentrators, Vivaspin 2, MWCO 10000 da	100 pz.
LLG09049762	Centrifugal concentrators, Vivaspin 2, MWCO 30000 da	25 pz.
LLG09049779	Centrifugal concentrators, Vivaspin 2, MWCO 30000 da	100 pz.
LLG09049763	Centrifugal concentrators, Vivaspin 2, MWCO 50000 da	25 pz.
LLG09049780	Centrifugal concentrators, Vivaspin 2, MWCO 50000 da	100 pz.
LLG09049764	Centrifugal concentrators, Vivaspin 2, MWCO 100000 da	25 pz.
LLG09049781	Centrifugal concentrators, Vivaspin 2, MWCO 100000 da	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

