



Condensers acc. to Allihn for extractors, DURAN tubing

Lenz

Condensers acc. to Allihn for extractors, DURAN® tubing

Made of DURAN® tubing. For condensing vapours. With cone, water connection with glass olives or GL 14 thread and two threaded plastic olives.



Code	Description	Packaging
LLG06204890	Condensers acc. to Allihn for extractors, DURAN® tubing, Extractor 30 ml, Effective length 160 mm, Cone 29/32 NS, Connection Glass olive, Cooling surface 100 cm ²	1 pz.
LLG06242774	Condensers acc. to Allihn for extractors, DURAN® tubing, Extractor 30 ml, Effective length 160 mm, Cone 29/32 NS, Connection GL 14, Cooling surface 100 cm ²	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

