



Extraction apparatuses acc. to Soxhlet, with Dimroth condenser, DURAN tubing

Lenz

Extraction apparatuses acc. to Soxhlet, with Dimroth condenser, DURAN® tubing

Extraction apparatus, Soxhlet-pattern. With Dimroth condenser and round-bottom flask.

Made of Duran® tubing.



Code	Description	Packaging
LLG09043001	Extraction apparatuses , Extractor 30 ml, Condenser 29/32 NS, Round bottom flasks Capacity 100 ml	1 pz.
LLG09043003	Extraction apparatuses , Extractor 100 ml, Condenser 45/40 NS, Round bottom flasks Capacity 250 ml	1 pz.
LLG09043005	Extraction apparatuses , Extractor 250 ml, Condenser 45/40 NS, Round bottom flasks Capacity 500 ml	1 pz.
LLG09043008	Extraction apparatuses , Extractor 200 ml, Condenser 45/40 NS, Round bottom flasks Capacity 250 ml	1 pz.
LLG09043004	Extraction apparatuses , Extractor 150 ml, Condenser 45/40 NS, Round bottom flasks Capacity 250 ml	1 pz.
LLG09043006	Extraction apparatuses , Extractor 500 ml, Condenser 60/46 NS, Round bottom flasks Capacity 1000 ml	1 pz.
LLG09043007	Extraction apparatuses , Extractor 1000 ml, Condenser 71/51 NS, Round bottom flasks Capacity 2000 ml	1 pz.
LLG09043009	Extraction apparatuses , Extractor 300 ml, Condenser 60/46 NS, Round bottom flasks Capacity 500 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

