

**Digital Abbe refractometer, DR-A1-Plus****ATAGO****Digital Abbe refractometer, DR-A1-Plus**

Very simple operation, where only the refraction boundary line has to be moved to the measurement cross hairs. The refractometer displays measurements directly (as a refractive index or in Brix (%)) in digits together with temperature. Measurements can be taken easily by anyone using this refractometer without taking an analogue scale reading in advance. Dispersive value can not be measured using this refractometer.

**Specifications**

**Measurement range:** 1.3000 - 1.7100 nD

**0.0 - 100.0 Brix**

**Resolution:** 0.0001 nD

**0.1 %Brix**

**Accuracy:** 0.0002 nD

**0.1 %Brix**

**Temperature range:** 0 to 50.0 °C

**Power supply:** a.c. adapter (100 to 240 V 50/60 Hz)

**Dimensions:** 130 x 290 x 310 mm

**Weight:** 6.0 kg



Code	Description	Packaging
LLG06283092	Digital Abbe refractometer, DR-A1-Plus, Type DR-A1-Plus	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

