

**Drum, Flask, STR4/2 for overhead shakers STR4****Stuart****Drum, Flask, STR4/2 for overhead shakers STR4**

Drum with two platforms for STR4. Accepts two 500ml conical flasks or reagent bottles. The support columns can be detached to give two platforms 10.5cm x 17.5cm each for separating funnels, jars or test tube racks.

With BioCote, silver-based, antimicrobial protection.

Drive Art. no 9.950 375 separately.



Code	Description	Packaging
LLG09950382	Drum, Flask, STR4/2 for overhead shakers STR4, Type STR4/2	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

