

**Drum figure-of-8 STR4/3 for overhead shakers STR4****Stuart****Drum figure-of-8 STR4/3 for overhead shakers STR4**

The drum has a platform which can be adjusted within the drum from horizontal to 45°. In this way the sample container can be rolled and turned end over end at the same time. Suitable for extractions, paint mixing, powders and granular samples. The sample container is held by an adjustable velcro strap. With BioCote, silver-based, antimicrobial protection.

Drive Art. no 9.950 375 separately.



Code	Description	Packaging
LLG09950384	Drum figure-of-8 STR4/3 for overhead shakers STR4, Type STR4/3	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

