


**Bottle holder STR4/4 for overhead shakers STR4****Stuart****Bottle holder STR4/4 for overhead shakers STR4**

Perspex four segment cradle with velcro straps for STR4. Holds 2 to 4 bottles each up to 1000ml capacity.

With BioCote, silver-based, antimicrobial protection.

Drive Art. no 9.950 375 separately.



Code	Description	Packaging
LLG09950385	Bottle holder STR4/4 for overhead shakers STR4, Type STR4/4	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

