



Wide-mouth reagent bottles, soda-lime glass

DWK Life Sciences

Wide-mouth reagent bottles, soda-lime glass

- To ISO 4796-2
- With ground joint neck
- With glass stopper
- Hydrolytic class 3



Code	Description	Packaging
LLG09071519	Wide-mouth reagent bottles, soda-lime glass, Colour amber, 500 ml, Diam. 86 mm, Height 164 mm, Neck 45/27 NS	1 pz
LLG09071408	Wide-mouth reagent bottles, soda-lime glass, Colour clear, 100 ml, Diam. 52 mm, Height 97 mm, Neck 29/22 NS	1 pz
LLG09071508	Wide-mouth reagent bottles, soda-lime glass, Colour amber, 100 ml, Diam. 52 mm, Height 97 mm, Neck 29/22 NS	1 pz
LLG09071424	Wide-mouth reagent bottles, soda-lime glass, Colour clear, 1000 ml, Diam. 107 mm, Height 200 mm, Neck 60/46 NS	1 pz
LLG09071524	Wide-mouth reagent bottles, soda-lime glass, Colour amber, 1000 ml, Diam. 107 mm, Height 200 mm, Neck 60/46 NS	1 pz
LLG09071414	Wide-mouth reagent bottles, soda-lime glass, Colour clear, 250 ml, Diam. 71 mm, Height 129 mm, Neck 34/24 NS	1 pz
LLG09071514	Wide-mouth reagent bottles, soda-lime glass, Colour amber, 250 ml, Diam. 71 mm, Height 129 mm, Neck 34/24 NS	1 pz
LLG09071405	Wide-mouth reagent bottles, soda-lime glass, Colour clear, 50 ml, Diam. 44 mm, Height 79 mm, Neck 24/20 NS	1 pz
LLG09071419	Wide-mouth reagent bottles, soda-lime glass, Colour clear, 500 ml, Diam. 86 mm, Height 164 mm, Neck 45/40 NS	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



