

**Magnetic stirrer with heating, 3-Position, SB162-3****Stuart****Magnetic stirrer with heating, 3-Position, SB162-3**

With 3 independently adjustable aluminium hotplates, Overheat protection mechanism and warning light which switches on when the plate temperature exceeds 50°C. With facility for mounting a 12.5mm retort rod at the rear of each plate and 3mm x 25mm PTFE coated followers.

With BioCote, silver-based, antimicrobial protection.

Specifications

Max. stirring quantity (H₂O): 15l

Speed range: 100 - 1500rpm

Heat output: 3 x 700W

Max. plate temperature: 325 °C

Set-up plate (W x D): 3 x 160 x 160 mm

Dimensions (W x D x H): 600 x 270 x 110mm

Weight: 11kg

Power supply: 230V, 50Hz

IP code: IP 32



Code	Description	Packaging
LLG09950424	Magnetic stirrer SB 162-3 with heating plate	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

