

**Laboratory bottles, DURAN®, with retrace code, with screw cap**

ISO 4796-1 compliant. Thread acc. to DIN 168-1 and approximate graduations. The ring mark on the shoulder of the 100ml to the 2000ml bottles makes the filling height of the nominal capacity visible. The 25ml bottle has a specially shaped glass rim for pouring (eliminating the need for an additional pouring ring). Neutral glass / Type I glass in compliance with USP/EP and JP. With easy to read scale and large labelling field for easy marking. Highly durable white ceramic. With Retrace code (Batch identification). Certificate available via the Internet. Autoclavable. With blue PP screw cap and pouring ring.



Code	Description	Packaging
LLG09072000	Laboratory bottles, DURAN®, with retrace code, with screw cap, 50 ml, Diam. 46 mm, Height with cap 91 mm, Neck thread 32 GL	1 pz
LLG09072005	Laboratory bottles, DURAN®, with retrace code, with screw cap, 500 ml, Diam. 86 mm, Height with cap 181 mm, Neck thread 45 GL	1 pz
LLG09072001	Laboratory bottles, DURAN®, with retrace code, with screw cap, 100 ml, Diam. 56 mm, Height with cap 105 mm, Neck thread 45 GL	1 pz
LLG09072010	Laboratory bottles, DURAN®, with retrace code, with screw cap, 1000 ml, Diam. 101 mm, Height with cap 230 mm, Neck thread 45 GL	1 pz
LLG09072004	Laboratory bottles, DURAN®, with retrace code, with screw cap, 10000 ml, Diam. 227 mm, Height with cap 415 mm, Neck thread 45 GL	1 pz
LLG09071680	Laboratory bottles, DURAN®, with retrace code, with screw cap, 150 ml, Diam. 62 mm, Height with cap 105 mm, Neck thread 45 GL	1 pz
LLG07076312	Laboratory bottles, DURAN®, with retrace code, with screw cap, 15000 ml, Diam. 268 mm, Height with cap 450 mm, Neck thread 45 GL	1 pz
LLG09072007	Laboratory bottles, DURAN®, with retrace code, with screw cap, 2000 ml, Diam. 136 mm, Height with cap 265 mm, Neck thread 45 GL	1 pz
LLG09072009	Laboratory bottles, DURAN®, with retrace code, with screw cap, 20000 ml, Diam. 288 mm, Height with cap 510 mm, Neck thread 45 GL	1 pz
LLG09072003	Laboratory bottles, DURAN®, with retrace code, with screw cap, 25* ml, Diam. 36 mm, Height with cap 74 mm, Neck thread 25 GL	1 pz
LLG09072002	Laboratory bottles, DURAN®, with retrace code, with screw cap, 250 ml, Diam. 70 mm, Height with cap 143 mm, Neck thread 45 GL	1 pz
LLG09071682	Laboratory bottles, DURAN®, with retrace code, with screw cap, 3500 ml, Diam. 160 mm, Height with cap 295 mm, Neck thread 45 GL	1 pz
LLG09072008	Laboratory bottles, DURAN®, with retrace code, with screw cap, 5000 ml, Diam. 182 mm, Height with cap 335 mm, Neck thread 45 GL	1 pz
LLG09071681	Laboratory bottles, DURAN®, with retrace code, with screw cap, 750 ml, Diam. 95 mm, Height with cap 203 mm, Neck thread 45 GL	1 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

