

**Tamper-evident caps for soda-lime square screw cap bottles, PP****Tamper-evident caps for soda-lime square screw cap bottles, PP**

Blue. For square screw cap bottles (soda-lime glass).

Type H: Narrow-neck, with high form thread

Type K: Wide-necked, with short form thread

Temperature (max.): 140°C



| Code        | Description                                                                                                             | Packaging |
|-------------|-------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG09072092 | Tamper-evident caps for soda-lime square screw cap bottles, PP, Type 32 K, Diam. 44 mm, Height 23 mm, Neck thread 32 GL | 1 pz.     |
| LLG09072041 | Tamper-evident caps for soda-lime square screw cap bottles, PP, Type 32 H, Diam. 45 mm, Height 32 mm, Neck thread 32 GL | 1 pz.     |
| LLG09072093 | Tamper-evident caps for soda-lime square screw cap bottles, PP, Type 45 K, Diam. 58 mm, Height 27 mm, Neck thread 45 GL | 1 pz      |
| LLG09072095 | Tamper-evident caps for soda-lime square screw cap bottles, PP, Type 60 K, Diam. 78 mm, Height 29 mm, Neck thread 60 GL | 1 pz      |
| LLG09072042 | Tamper-evident caps for soda-lime square screw cap bottles, PP, Type 45 H, Diam. 60 mm, Height 35 mm, Neck thread 45 GL | 1 pz      |
| LLG09072094 | Tamper-evident caps for soda-lime square screw cap bottles, PP, Type 54 K, Diam. 69 mm, Height 29 mm, Neck thread 54 GL | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

