

**Quantitative filter paper, grade 40****Whatman****Quantitative filter paper, grade 40**

Circles. Quantitative. Ashless.

Flow rate:

Retention:

Wet strength:

Thickness:

Ash: Medium

Medium

Low

0.21mm

0.007%



| Code        | Description                                        | Packaging |
|-------------|----------------------------------------------------|-----------|
| LLG09056403 | Quantitative filter paper, grade 40, Diam. 110 mm  | 100 pz.   |
| LLG09056404 | Quantitative filter paper, grade 40, Diam. 125 mm  | 100 pz.   |
| LLG09056405 | Quantitative filter paper, grade 40, Diam. 150 mm  | 100 pz.   |
| LLG09056406 | Quantitative filter paper, grade 40, Diam. 185 mm  | 100 pz.   |
| LLG09056407 | Quantitative filter paper, grade 40, Diam. 240 mm  | 100 pz.   |
| LLG09049035 | Quantitative filter paper, grade 40, Diam. 32 mm   | 100 pz.   |
| LLG09056399 | Quantitative filter paper, grade 40, Diam. 42.5 mm | 100 pz.   |
| LLG09056429 | Quantitative filter paper, grade 40, Diam. 47 mm   | 100 pz.   |
| LLG09056400 | Quantitative filter paper, grade 40, Diam. 55 mm   | 100 pz.   |
| LLG09056401 | Quantitative filter paper, grade 40, Diam. 70 mm   | 100 pz.   |
| LLG09056402 | Quantitative filter paper, grade 40, Diam. 90 mm   | 100 pz.   |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

