

**Quantitative filter paper, grade 41****Whatman****Quantitative filter paper, grade 41**

Circles. Quantitative. Ashless. For filtration of coarse and gelatinous precipitates.

Supplied in box of 100.

Flow rate:

Retention:

Wet strength:

Thickness:

Ash: Fast

Medium

Low

0.22mm

0.007%



Code	Description	Packaging
LLG09056413	Quantitative filter paper, grade 41, Diam. 110 mm	100 pz.
LLG09056414	Quantitative filter paper, grade 41, Diam. 125 mm	100 pz.
LLG09056415	Quantitative filter paper, grade 41, Diam. 150 mm	100 pz.
LLG09056416	Quantitative filter paper, grade 41, Diam. 185 mm	100 pz.
LLG09056417	Quantitative filter paper, grade 41, Diam. 240 mm	100 pz.
LLG09056410	Quantitative filter paper, grade 41, Diam. 55 mm	100 pz.
LLG09049045	Quantitative filter paper, grade 41, Diam. 60 mm	100 pz.
LLG09056411	Quantitative filter paper, grade 41, Diam. 70 mm	100 pz.
LLG09056412	Quantitative filter paper, grade 41, Diam. 90 mm	100 pz.
LLG09056563	Quantitative filter paper, grade 41, Diam. mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

