



Hydrometers, relative density (S.G.)

Hydrometers, relative density (S.G.)

Shot poised. Without thermometer. Reference temperature: +20°C.



Code	Description	Packaging
LLG09004008	Hydrometers, relative density (S.G.), Length 280 mm	1 pz.
LLG09004012	Hydrometers, relative density (S.G.), Length 280 mm	1 pz.
LLG09004014	Hydrometers, relative density (S.G.), Length 280 mm	1 pz.
LLG09004018	Hydrometers, relative density (S.G.), Length 300 mm	1 pz.
LLG09004057	Hydrometers, relative density (S.G.), Length 300 mm	1 pz.
LLG09004060	Hydrometers, relative density (S.G.), Length 300 mm	1 pz.
LLG09004065	Hydrometers, relative density (S.G.), Length 300 mm	1 pz.
LLG09004006	Hydrometers, relative density (S.G.), Length 280 mm	1 pz.
LLG09004010	Hydrometers, relative density (S.G.), Length 280 mm	1 pz.
LLG09004016	Hydrometers, relative density (S.G.), Length 280 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

