



Membrane filters, Nylon

Membrane filters, Nylon

Hydrophilic, durable and compatible with a wide range of aqueous solutions, organic solvents and for the filtration of biological preparations and media. Autoclavable at 121 °C.



Code	Description	Packaging
LLG09951042	Membrane filters, Nylon, Diam. 13 mm, Pore size 0.2 µm	100 pz.
LLG09951045	Membrane filters, Nylon, Diam. 13 mm, Pore size 0.45 µm	100 pz.
LLG09951046	Membrane filters, Nylon, Diam. 25 mm, Pore size 0.45 µm	100 pz.
LLG09951043	Membrane filters, Nylon, Diam. 47 mm, Pore size 0.2 µm	100 pz.
LLG09951047	Membrane filters, Nylon, Diam. 47 mm, Pore size 0.45 µm	100 pz.
LLG09951049	Membrane filters, Nylon, Diam. 47 mm, Pore size 0.8 µm	100 pz.
LLG09951044	Membrane filters, Nylon, Diam. 90 mm, Pore size 0.2 µm	50 pz.
LLG09951048	Membrane filters, Nylon, Diam. 90 mm, Pore size 0.45 µm	50 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

