



Disposable Filtration Capsules, Polycap HD, PP

Disposable Filtration Capsules, Polycap HDä, PP

PP monofilament depth filters for use as prefilters or for the filtration of samples in water or organic solvents. The high loading capacity allows the purification of solutions with a high particle load.

- Choice of filtration areas and porosities
- Maximum operating pressure 4.1 bar forward, 3.0 bar reverse
- Capsules may be used in either direction
- Low hold-up volume with air purge
- Autoclavable
- Operating temperature max. 60 °C
- Connectors: stepped barb (SB), hose barb (HB) and female national pipe thread (FNPT)
- Filter PP



Code	Description	Packaging
LLG09951118	Disposable Filtration Capsules, Polycap HD, PP, Porosity 1 µm, Filter diam. 36 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz.
LLG09951122	Disposable Filtration Capsules, Polycap HD, PP, Porosity 1 µm, Filter diam. 75 mm, Connection in / out 3/8-1/2" SB / 3/8-1/2" SB, Sterile -	1 pz.
LLG09951125	Disposable Filtration Capsules, Polycap HD, PP, Porosity 20 µm, Filter diam. 75 mm, Connection in / out 3/8-1/2" SB / 3/8-1/2" SB, Sterile -	1 pz.
LLG09951120	Disposable Filtration Capsules, Polycap HD, PP, Porosity 10 µm, Filter diam. 36 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz.
LLG09951121	Disposable Filtration Capsules, Polycap HD, PP, Porosity 20 µm, Filter diam. 36 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz.
LLG09951119	Disposable Filtration Capsules, Polycap HD, PP, Porosity 5 µm, Filter diam. 36 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz.
LLG09951123	Disposable Filtration Capsules, Polycap HD, PP, Porosity 5 µm, Filter diam. 75 mm, Connection in / out 3/8-1/2" SB / 3/8-1/2" SB, Sterile -	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

