



## Disposable Filtration Capsules, Polycap TF

### Disposable Filtration Capsules, Polycap TFä

Hydrophobic PTFE membrane for applications involving organic solvents and chemically aggressive solutions.

- Maximum operating pressure of 4.1 bar forward and 3.0 bar reverse
- 1 µm capsules incorporate a PP prefilter
- Autoclavable up to 132 °C maximum
- Connectors: stepped barb (SB) and female national pipe thread (FNPT)



Code	Description	Packaging
LLG09951126	Disposable Filtration Capsules, Polycap TF, Porosity 0.2 µm, Filter diam. 36 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz
LLG09951127	Disposable Filtration Capsules, Polycap TF, Porosity 1.0 µm, Filter diam. 36 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz
LLG09951128	Disposable Filtration Capsules, Polycap TF, Porosity 0.2 µm, Filter diam. 75 mm, Connection in / out 1/4-3/8" SB / 1/4-3/8" SB, Sterile -	1 pz
LLG09951129	Disposable Filtration Capsules, Polycap TF, Porosity 1.0 µm, Filter diam. 75 mm, Connection in / out 3/8-1/2" SB / 3/8-1/2" SB, Sterile -	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

