


**Light Boxes slimlite plano****Kaiser Fototechnik****Light Boxes slimlite plano**

Super flat light boxes, only 8 mm thin, with built-in rechargeable battery. Illumination by daylight-type LEDs (approx. 5000 Kelvin). Touch switch for ON/OFF and brightness control, cm and inch scales. Optional mains or battery operation. Lithium ion battery: 1050 mA. Mains/charging adapter included.



Code	Description	Packaging
LLG06270828	Light Boxes slimlite plano, Illuminated area (W x D) 220 x 160 mm, Lighting 23 LEDs	1 pz.
LLG06270829	Light Boxes slimlite plano, Illuminated area (W x D) 320 x 228 mm, Lighting 33 LEDs	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

