


**Digital block thermostats QB, for QBD, QBH series****Grant****Digital block thermostats QB, for QBD, QBH series**

Accessory blocks for QB series only. Interchangeable. Aluminium. Block dimensions (L x W x H): 140 x 50 x 63 mm. For e.g. test tubes or microliter vials from 0.2 ml to 2.0 ml.



Code	Description	Packaging
LLG09951666	Digital block thermostats QB, for QBD, QBH series, Type QB-H, For tubes diam. 0.2 ml mm, No.of holes 56	1 pz
LLG09951665	Digital block thermostats QB, for QBD, QBH series, Type QB-24, For tubes diam. 24 mm mm, No.of holes 5	1 pz
LLG09951667	Digital block thermostats QB, for QBD, QBH series, Type QB-E0, For tubes diam. 0.5 ml mm, No.of holes 24	1 pz
LLG09951668	Digital block thermostats QB, for QBD, QBH series, Type QB-E1, For tubes diam. 1.5 ml mm, No.of holes 24	1 pz
LLG09951661	Digital block thermostats QB, for QBD, QBH series, Type QB-10, For tubes diam. 10 mm mm, No.of holes 24	1 pz
LLG09951663	Digital block thermostats QB, for QBD, QBH series, Type QB-16, For tubes diam. 16 mm mm, No.of holes 12	1 pz
LLG09951664	Digital block thermostats QB, for QBD, QBH series, Type QB-18, For tubes diam. 18 mm mm, No.of holes 12	1 pz
LLG09951669	Digital block thermostats QB, for QBD, QBH series, Type QB-E2, For tubes diam. 2.0 ml mm, No.of holes 24	1 pz
LLG06266809	Digital block thermostats QB, for QBD, QBH series, Type QB-E5, For tubes diam. 5.0 ml mm, No.of holes 12	1 pz.
LLG06075814	Digital block thermostats QB, for QBD, QBH series, Type QB-50, For tubes diam. 50 ml mm, No.of holes 4	1 pz.
LLG09951670	Digital block thermostats QB, for QBD, QBH series, Type QB-0, For tubes diam. solid, for user to drill as required mm, No.of holes	1 pz
LLG09951662	Digital block thermostats QB, for QBD, QBH series, Type QB-12, For tubes diam. 12 mm mm, No.of holes 24	1 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

