




Test Sieves, 200 x 50mm

Test Sieves, 200 x 50mm

Features:

- Previously unattained product quality owing to advanced manufacturing processes and continuous optical inspection.
 - High quality joining technology guarantees permanently taut wire mesh.
 - High level of corrosion resistance and easy to clean owing to the use of high-alloy stainless steel.
 - "one-piece cast" sieve frame and seamless mesh joint prevent cross contamination.
 - Maximum stability and optimum density of the sieve stack.
 - 15% lighter than traditional sieves whilst at the same time having larger sieving area.
 - Clear sieve marking and traceability by individual laser engraving on each sieve.
- DIN ISO 3310-1. Dimensions (dia. x h): 200mm x 50mm. **Test sieves with height 25mm available on request.**



Code	Description	Packaging
LLG09226692	Test sieves 200x50mesh size 1,00mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

