

**Test Sieves, 200 x 25 mm****Retsch****Test Sieves, 200 x 25 mm****Features:**

- Previously unattained product quality owing to advanced manufacturing processes and continuous optical inspection.
 - High quality joining technology guarantees permanently taut wire mesh.
 - High level of corrosion resistance and easy to clean owing to the use of high-alloy stainless steel.
 - "one-piece cast" sieve frame and seamless mesh joint prevent cross contamination.
 - Maximum stability and optimum density of the sieve stack.
 - 15% lighter than traditional sieves whilst at the same time having larger sieving area.
 - Clear sieve marking and traceability by individual laser engraving on each sieve.
- DIN ISO 3310-1. Dimensions (dia. x H): 200mm x 25mm.



Code	Description	Packaging
LLG09226765	Test Sieves, 200 x 25 mm, Mesh size 100 µm	1 pz.
LLG09226766	Test Sieves, 200 x 25 mm, Mesh size 106 µm	1 pz.
LLG09226767	Test Sieves, 200 x 25 mm, Mesh size 112 µm	1 pz.
LLG09226768	Test Sieves, 200 x 25 mm, Mesh size 125 µm	1 pz.
LLG09226769	Test Sieves, 200 x 25 mm, Mesh size 140 µm	1 pz.
LLG09226770	Test Sieves, 200 x 25 mm, Mesh size 150 µm	1 pz.
LLG09226771	Test Sieves, 200 x 25 mm, Mesh size 160 µm	1 pz.
LLG09226772	Test Sieves, 200 x 25 mm, Mesh size 180 µm	1 pz.
LLG09226750	Test Sieves, 200 x 25 mm, Mesh size 20 µm	1 pz.
LLG09226773	Test Sieves, 200 x 25 mm, Mesh size 200 µm	1 pz.
LLG09226774	Test Sieves, 200 x 25 mm, Mesh size 212 µm	1 pz.
LLG09226775	Test Sieves, 200 x 25 mm, Mesh size 224 µm	1 pz.
LLG09226751	Test Sieves, 200 x 25 mm, Mesh size 25 µm	1 pz.
LLG09226776	Test Sieves, 200 x 25 mm, Mesh size 250 µm	1 pz.
LLG09226777	Test Sieves, 200 x 25 mm, Mesh size 280 µm	1 pz.
LLG09226778	Test Sieves, 200 x 25 mm, Mesh size 300 µm	1 pz.
LLG09226779	Test Sieves, 200 x 25 mm, Mesh size 315 µm	1 pz.
LLG09226752	Test Sieves, 200 x 25 mm, Mesh size 32 µm	1 pz.
LLG09226780	Test Sieves, 200 x 25 mm, Mesh size 355 µm	1 pz.
LLG09226753	Test Sieves, 200 x 25 mm, Mesh size 36 µm	1 pz.



Code	Description	Packaging
LLG09226754	Test Sieves, 200 x 25 mm, Mesh size 38 µm	1 pz.
LLG09226755	Test Sieves, 200 x 25 mm, Mesh size 40 µm	1 pz.
LLG09226781	Test Sieves, 200 x 25 mm, Mesh size 400 µm	1 pz.
LLG09226782	Test Sieves, 200 x 25 mm, Mesh size 425 µm	1 pz.
LLG09226756	Test Sieves, 200 x 25 mm, Mesh size 45 µm	1 pz.
LLG09226783	Test Sieves, 200 x 25 mm, Mesh size 450 µm	1 pz.
LLG09226757	Test Sieves, 200 x 25 mm, Mesh size 50 µm	1 pz.
LLG09226784	Test Sieves, 200 x 25 mm, Mesh size 500 µm	1 pz.
LLG09226758	Test Sieves, 200 x 25 mm, Mesh size 53 µm	1 pz.
LLG09226759	Test Sieves, 200 x 25 mm, Mesh size 56 µm	1 pz.
LLG09226785	Test Sieves, 200 x 25 mm, Mesh size 560 µm	1 pz.
LLG09226786	Test Sieves, 200 x 25 mm, Mesh size 600 µm	1 pz.
LLG09226760	Test Sieves, 200 x 25 mm, Mesh size 63 µm	1 pz.
LLG09226787	Test Sieves, 200 x 25 mm, Mesh size 630 µm	1 pz.
LLG09226761	Test Sieves, 200 x 25 mm, Mesh size 71 µm	1 pz.
LLG09226788	Test Sieves, 200 x 25 mm, Mesh size 710 µm	1 pz.
LLG09226762	Test Sieves, 200 x 25 mm, Mesh size 75 µm	1 pz.
LLG09226763	Test Sieves, 200 x 25 mm, Mesh size 80 µm	1 pz.
LLG09226789	Test Sieves, 200 x 25 mm, Mesh size 800 µm	1 pz.
LLG09226790	Test Sieves, 200 x 25 mm, Mesh size 850 µm	1 pz.



Code	Description	Packaging
LLG09226764	Test Sieves, 200 x 25 mm, Mesh size 90 µm	1 pz.
LLG09226791	Test Sieves, 200 x 25 mm, Mesh size 900 µm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

